



# Low Beds in Healthcare: Proven Success, but When Is Low Too Low?

**Balancing Fall Prevention, Patient Safety, and Operational Efficiency**

---

**Abstract:** Low beds are widely recognized for their role in reducing fall-related injuries within healthcare environments. However, determining the optimal bed height is essential—lowering beds too much can introduce new risks impacting patient safety, caregiver ergonomics, and infection control (Merryweather AS, 2015). This paper synthesizes clinical evidence, biomechanical principles, and operational realities to guide decision-makers in finding the right balance between fall prevention, egress safety, and workflow efficiency.



## THE IDEAL TRANSFER HEIGHT

Research consistently identifies 18 inches as the ideal average transfer height, aligning with the popliteal height for most patients. This typically corresponds to a 12-inch bed frame combined with a standard 6-inch foam surface. Speed in adjusting the bed to this optimal height is critical: delays can increase the likelihood of unsafe transfers or lead caregivers to omit necessary adjustments, thereby compromising patient safety (Christman M, 2015; Merryweather AS, 2015).



\*\*It is recommended that all residents be assessed for their ideal egress height\*\*


## THE CASE FOR LOW BEDS

Lowering bed height is an effective strategy for reducing the severity of injuries from rollouts, particularly head and pelvic injuries. Biomechanical analysis reveals that impact velocity, and thus injury risk, decreases as bed height is reduced; however, the benefit plateaus below a certain point (Bowers B, 2008).



**Evidence shows:**  
**Lowering bed height**  
reduces **injury risk**,  
but **benefits plateau**  
below **7 inches**.

### Reduction in Impact by Frame Configuration Change



| BED CONFIGURATION | IMPACT (Ns) |
|-------------------|-------------|
| 10.1" Frame       | 387         |
| 7" Frame          | 122         |
| 3.9" Frame        | 39          |

*The largest reduction in impact occurs between the 10.1" and 7" frame configurations. Subsequent reductions are smaller, indicating diminishing returns at lower frame heights.*



## WHEN LOWER BECOMES A LIABILITY

There is a critical point below which further lowering the bed system yields diminishing returns. Ultra-low beds introduce several operational and safety risks:



**Infection Control:**  
Beds closer to the floor increase contamination risk.



**Egress Challenges:**  
Patients find it harder to exit ultra-low beds, heightening the risk of falls during transfers.



**Caregiver Ergonomics:**  
Ultra-low beds are difficult to raise to a safe working height, increasing the risk of caregiver injury due to frequent bending.



**Operational Delays:**  
Waiting for beds to reach ideal egress height can slow down care delivery.

As highlighted in the supporting context, the speed of reaching the ideal egress height is just as critical as the minimum height itself. Beds that are too low can complicate both patient movement and caregiver workflow, sometimes negating the intended safety benefits (Morse JM, 2015).

## BEYOND HEIGHT: HOLISTIC BED SYSTEM DESIGN

Optimal bed safety encompasses multiple factors beyond height reduction. While lowering bed height can mitigate injury severity, the protective effect plateaus below approximately 7 inches. Accurate definitions of falls, including rollouts, are essential for valid reporting and targeted interventions.

Excessive lowering may introduce secondary risks, such as increased contamination potential, compromised caregiver ergonomics, and delayed care delivery. Furthermore, premium “floor-level” designs frequently lack demonstrable outcome improvements in clinical practice.





## CONCLUSION

Low beds are an important tool in fall prevention strategies, but lowering bed height beyond a certain point offers negligible benefits and introduces new risks. Evidence shows:

- Going below 7 inches is less effective than alternative solutions when calculating impact velocity.
- Ultra-low designs can increase infection risk, neglect ergonomic, and egress risks.
- Premium “floor-level” features often do not deliver meaningful outcome improvements.

**RECOMMENDATION:** Select bed systems that strike a balance between low height and rapid adjustment to optimal egress height, incorporate ergonomic design for caregivers, and feature high-quality surfaces. Prioritize evidence-based features over marketing claims and ensure a holistic approach to bed safety and patient outcomes.

## WORKS CITED

Becker, C. (2010). Fall prevention in nursing homes. *Clinics in Geriatric Medicine*, 26(4), 693–704.

Boushon, B., & Graber, N. (2012). How-to guide: Reducing patient injuries from falls. *Institute for Healthcare Improvement*.

Bowers, B., & Johnson-Campbell, L. (2008). Biomechanical evaluation of injury severity associated with patient falls from bed. *Rehabilitation Nursing*, 33(6), 253–259.

Christman, M., & Miller, J. (2015). Analysis of the influence of hospital bed height on kinematic parameters associated with patient falls during egress. *Procedia Manufacturing*, 3, 280–287.

Crane, B. C. (2011). A floor covering to reduce injury in hospital falls: A multidisciplinary approach. Platform presentation at the Connecticut Chapter Physical Therapy Association Annual Conference.

Merryweather, A. S., & Miller, J. (2015). Effects of bed height on the biomechanics of hospital bed entry and egress. *Work*, 52(3), 707–713.

Morse, J. M., & Evans, A. (2015). The safety of hospital beds: Ingress, egress, and in-bed mobility. *Global Qualitative Nursing Research*, 2, 1–20.

Raymond, D. E., & Evans, A. (2011). Biomechanics and injury risk assessment of falls onto protective floor mats. *Rehabilitation Nursing*, 36(6), 248–254.

## ABOUT JOERNS HEALTHCARE

Joerns Healthcare provides healthcare technology and related services for patient care across the post-acute continuum. It is dedicated to simplifying this process by providing practitioners with the products, technology and data they need, allowing them to concentrate on patient care. Joerns also offers over 130 years of experience and a national service footprint that help clinicians navigate the complexities of modern accountable care while optimizing patient outcomes and cost.

[info@joerns.com](mailto:info@joerns.com)

800.826.0270