

Joerns® Bari10A Bed



leading bed evolution

# Joerns® Bari10A Bed Leading Bariatric Solutions

The new Joerns Bari10A bed was specifically designed to address bariatric patient dignity while providing patient safety and comfort when patient weight exceeds standard surface weight limits. The Bari10A bariatric bed is designed with a fully articulated mattress support platform with integrated width extension for comfort and transport modes. The Bari10A bed design allows for independent head and knee movement, up/down adjustments, gravity assist re-positioning, while ensuring optimal clinical outcomes for your bariatric patients.

## Classifications:

- UL60601-1, CAN/CSA C22.2 NO. 601.1, IEC 60601-2-38, CAN/CSA C22.2 NO. 601.2.38

## Joerns® Bari10A bed features:

- 1 The safe working load is 1050 lbs.
- 2 Integrated battery backup system in the event of power outage or patient transfer.
- 3 Easy to operate integrated scale system for resident weight monitoring.
- 4 Pendant control for all bed functions with separate selective lock-out to override functions.
- 5 36" wide mattress support platform sections easily expands to 48" width without tools while using one hand.
- 6 80" mattress support platform, with an integrated foot end extension section extending to an extra 8" for taller patients.
- 7 Easily configurable in four sizes – 36" X 80", 48" X 80", 48" X 88", 36" x 88".
- 8 18<sup>3/4</sup>" to 33<sup>3/4</sup>" travel range offers flexibility for any care need, allowing for patient and caregiver safety and ergonomic benefits.
- 9 Attractive, non-bariatric style bed panels in a soft tone frame color.
- 10 Side rails on head and foot sections fold down for easy patient access and caregiver safety.
- 11 Four DC motors provide quiet and smooth operation.
- 12 Heavy-duty casters for easier mobility, featuring centerline caster locks.
- 13 Four-section sleep surface is infinitely positionable while helping with the reduction of skin shear.
- 14 No tools are needed to remove the head and foot panels.
- 15 Optional Fifth Wheel enables safe and efficient one person transport.



Joerns Healthcare understands the healthcare bed market. As a leading manufacturer of healthcare beds, we have developed the Bari10A bariatric bed to complete your healthcare bed needs. The new Joerns Bari10A bed is designed for a 1050 lb. work load and is sturdy throughout, with four motors providing independent head and knee, up/down and Trendelenburg and reverse Trendelenburg functions. We are confident in the design, the production and strength of our new Bari10A bed and guarantee the Joerns Bari10A bed for a period of one year from date of delivery.

## Trendelenburg/ Reverse Trendelenburg



Trendelenburg and Reverse Trendelenburg positioning allows staff to provide additional angles of care for specific patient needs.

## Safe Floor Gap



The lowered side rails maintain a safe floor gap to keep staff safe and avoid potential foot injury when lowering the bed.

## Low/High Position



Fully articulated mattress support platform can be used in the low and high positions to ensure safety while allowing for safe patient handling by the caregiver.

## Extended Width

Easy, convenient, one-hand width extensions provide for restraint-free support for patient transfers.



## Extended Length

Integrated width and length extension section makes transports more seamless.



## 5th Wheel



The optional Fifth Wheel assists in mobility to facilitate safe and efficient one-person patient transport.

## Scale Interface



Integrated electronic bed scale system, provides accurate, repeatable weights to be accessed by staff.

## Patient Helper



Optional Patient Helper has an adjustable front and back that assists in patient repositioning and ingress and egress. Designed specifically for bariatric patients.

# Joerns® Bari10A Bed Leading Bariatric Solutions

## Side Rail Release



Side rails on head and foot sections easily fold down for quick patient access and caregiver safety.

## CPR Release



Emergency manual release for head section.

## Angle Indicator



The angle indicator allows caregivers to read the TRT angle.

## Lockout Function

Separate selective override lock-out function for caregiver to limit bed position for specific patient care.



## Joerns Therapeutic Foam and Low Air Loss Mattress Solutions

### Powered and Non-Powered Bariatric Support Surfaces



The Bari10A features two optional support surfaces. Both surfaces are designed to provide comfort and protection for the obese or morbidly obese patients weighing up to 1000 lbs. and at risk of developing pressure ulcers.

Based on a thorough knowledge of our customers' needs, Joerns' therapeutic support mattresses were designed to make the caregiver's job easier and safer, while ensuring optimal clinical results.

The BariMatt 1000 pressure redistribution surface is a bariatric foam mattress which is specifically designed for the Bari10A frame while providing patients with a safe, comfortable sleep surface.

The Arise 1000 is a true low air loss solution that wicks away moisture and prevents skin maceration. With three modes of operation and eight therapeutic comfort control settings, the Arise 1000 is ideal for wound prevention and treatment of pressure ulcers up to stage IV.

Joerns Healthcare is committed to continuing our long tradition of caring for your staff and your patients.



Joerns Healthcare Inc.  
5001 Joerns Drive, Stevens Point WI 54481  
(P) 800.826.0270 • (F) 800.457.8827  
www.joerns.com • email: info@joerns.com