

	Max Weight Capacity	Lift Range of Boom*	Width	Length	Base Clearance	Finish	Shipping Weight
<b>Manual Hydraulic</b>							
HML400	400 lbs. (182 kg)	28"-77" (71-196 cm)	24"-40.5" (61-103 cm)	43.5" (110 cm)	5.5" (14 cm)	Painted Tan	96 lbs. (44 kg)
C-HLA-2	400 lbs. (182 kg)	29"-74.5" (74-189 cm)	24"-39.5" (61-100 cm)	41" (104 cm)	5.5" (14 cm)	Chrome	95 lbs. (43 kg)
C-HLA-2T	400 lbs. (182 kg)	37"-82.5" (94-210 cm)	24"-39.5" (61-100 cm)	41" (104 cm)	5.5" (14 cm)	Chrome	95 lbs. (43 kg)
C-CB-L2	400 lbs. (182 kg)	28.5"-74" (72-188 cm)	23" (58 cm)	39" (99 cm)	4.5" (11 cm)	Chrome	95 lbs. (43 kg)
C-HLA	400 lbs. (182 kg)	30.5"-78" (77-198 cm)	24"-39.5" (61-100 cm)	41" (104 cm)	5.5" (14 cm)	Chrome	84 lbs. (38 kg)
HOY-ADVANCE-H	340 lbs. (155 kg)	21.5"-73.5" (55-187 cm)	25"-44.5" (64-113 cm)	50" (127 cm)	4.5" (11 cm)	Painted Grey	63.6 lbs. (28.9 kg)
<b>Power Lifters</b>							
HPL402	400 lbs. (182 kg)	28"-77" (71-196 cm)	24"-40.5" (61-103 cm)	43.5" (110 cm)	5" (13 cm)	Painted Tan	114 lbs. (52 kg)
P-C-HLA-2	400 lbs. (182 kg)	29"-74.5" (74-189 cm)	24"-39.5" (61-100 cm)	41" (104 cm)	5.5" (14 cm)	Chrome	112 lbs. (51 kg)
HPL600WB	600 lbs. (273 kg)	23"-63" (58-160 cm)	26"-47" (66-119 cm)	49.375" (125 cm)	4.5" (10 cm)	Painted Tan	177 lbs. (80 kg)
HOY-ADVANCE-E	340 lbs. (155 kg)	21.5"-73.5" (55-187 cm)	25"-44.5" (64-113 cm)	50" (127 cm)	4.5" (11 cm)	Painted Grey	69.7 lbs. (31.7 kg)
PRO300	297 lbs. (135 kg)	44"-64" (112-163 cm)	24"-48" (61-122 cm)	46"-48" (117-122 cm)	4.0" (10 cm)	Painted White	132 lbs. (60 kg)
*Powered stand-aid includes (1) medium and (1) large long back sling, (1) battery, (1) charger							
<b>Pool Lifter</b>							
SS-HSP-1	400 lbs. (182 kg)	62" (157 cm)				Stainless Steel	66 lbs. (30 kg)
*Floor to pivot point of cradle attachment.							

### Hoyer Digital Scales

- Fast, easy to use and accurate
- Simple push-button operation
- Digital readout up to 600 lbs. (273 kg)
- Battery life: Approximately 3,000 readings
- Powered by 9v battery

### Model 59011A fits:

- Manual: HML400, C-HLA-2
- Power: HPL402, P-C-HLA-2

### Model HLS-5-D fits:

- C-HLA

### Model HPL600SC fits:

- HPL600WB

Note: 600 lb. capacity of scale will not increase lifting ability of 400 lb. capacity lifters.