



Hoyer Single Patient Slings

- Made of soft Polypropylene, the same material used by physicians for surgical gowns
- Patient Specific – slings are designed for single patient use and discarded when sling is no longer required
- Hygienic - Reduces risk of cross contamination between shared slings
- Eliminates need to launder patient slings
- Color-coded trim for easy size identification
- Color-coded loops for consistent attachment to cradle/lift
- Weight capacity of 500 lbs.
- Available as (1) each or in cartons of 10 slings
- You must use Hoyer slings with Hoyer lifts

Available in 2 configurations:

1. 6-Point Sling Configuration: Medium and Large Full Back Slings.
 2. 4-Point Sling Configuration: Medium and Large Comfort Slings.
- Hoyer Single Patient Slings cannot be washed. Once the sling requires laundering, you must dispose of sling.

NA25520	Comfort	Medium, each
NA25521	Comfort	Large, each
NA25523	Full Back	Medium, each
NA25524	Full Back	Large, each



* Warranty covers foam for material defects and workmanship.
Cover warranties vary. Call 800.826.0270 for details.



Hoyer Slide Sheets

Made of ultra low friction material

- Washable
- Available in three sizes for different transfer applications.
 1. Small sheets are to be used under resident's feet along with a large sheet when moving the resident in prone position.
 2. Medium sheets are to be used for positioning the resident's hips and mid section. Medium sheets can also be used for repositioning on smaller/ shorter residents in the prone position.
 3. Large sheets are to be used for lateral repositioning of average to tall residents in a prone position.



Small & medium sheets shown above.

Hoyer Slide Sheets do not have handles.

1. Handles promote poor ergonomic posture by forcing caregiver to grip the slide sheet by the handle.
2. Handles have the potential of harboring infectious disease.
3. High temperature washing and drying can cause handles to distort or shrink, making the handles difficult to use.

NA25525	Small	24" x 18"
NA25526	Medium	36" x 25"
NA25527	Large	60" x 31"

