




Both resident and caregiver benefit from the new Deluxe Assist Handle

Enhance Resident Safety with NEW Deluxe Assist Handle

- New shrouded design minimizes limb entrapment and entanglement while featuring two locking positions for additional safety
- Meets the current IEC 60601-2-38 requirements in all positions. Compliant with dimensions for zones 1, 2, 3 and 4 of the March 2006 FDA Entrapment Guidelines. Also addresses areas of concern identified as zones 5, 6, and 7 in the March 2006 FDA Entrapment Guidelines
- Mounted at an ergonomic height to best assist in safe ingress/egress while promoting independent repositioning
- Features ergonomic and sturdy yet soft touch handholds for improved resident grip and comfort
- Easy installation and quick removal without any tools needed, minimizes time demands on staff
- Optional integrated pendant holder on both sides making pendant stowage and retrieval simple for both residents and caregivers
- For use with EasyCare® 5, 7, 9, UltraCare® and UltraCare XT bed models



Industry's best assist handle for resident egress/ingress promoting safety and comfort. Shown in two locking positions.



Deluxe assist handle stays with the bed frame and resident.



Easily removed for storage.



Deluxe assist handle shown with both pendant holder styles.



Enhance safety and comfort for your residents and caregivers while promoting independence.

Order the Deluxe Assist Handle today! Call 800.826.0270.



**Joerns
Healthcare**

A Tradition of Caring • www.joerns.com

6180078 RevB • 09-0337