

BioClinic Surgery, Recovery and General  
Foam Positioners



transforming wound management

## Clinically-sound, cost-effective solutions for the management of diverse patient positioning needs and protection against peripheral nerve and skin injury

In today's hospital environment, positioners are needed for patient protection and comfort and are integral to many operating room procedures. The BioClinic line of foam positioners and protectors offers foam density options and strategic product designs to match patient requirements in most operating room and care setting environments. BioClinic foam positioners provide cost effective clinical solutions to improve pressure redistribution and enhance the care provided in hospitals, nursing homes and home care settings.

### Choose from:

- **Standard 1.2 lb. density**  
Provides excellent economical pressure redistribution and support
- **Premium 1.5 lb. density**  
Offers superior protection and pressure redistribution

### Features and Benefits:

- Foam density
- Highly resilient foam provides excellent patient support
- Non-toxic, latex free
- Hypoallergenic and non-irritating
- Flame-retardant
- Meets CAL 117
- Compression packaging
- Many products are vacuum-packed to reduce facility storage space needs
- Single patient use reduces risk of cross-contamination

### Head Cradles

Positioner provides comfortable immobilization of head and neck. May be used for lateral, supine or prone positions. Available in four sizes. Adult size also available in viscoelastic-topped frame design providing additional comfort and support.

Foam	Adult	Adult Slotted	Adult Large Frame	Adult Large Frame Visco Slotted	Child	Pediatric
● <b>Standard</b>	8815	n/a	n/a	n/a	8814	8813
● <b>Premium</b>	9817	9824	9818	550269	9814	9813



### Round Foam Ring Positioners

A versatile positioner that provides anatomical support and comfort. Available in a variety of sizes and foam densities.

Foam	5"	7"	9"	All-in-One
● <b>Standard</b>	8804	8807	8809	n/a
● <b>Premium</b>	8805	n/a	550211-CC	9809



## Ulnar Nerve Protectors

Protects against pinching and entrapment of the ulnar nerve. Provides cushioning on bony prominences of the arm. One size fits most, to reduce inventory needs.

- Standard 4816
- Premium 9816



## Heel/Ankle Protectors

Comfortable and easy to use. Provides cushioning of the foot, ankle and heel areas.

- Standard 4815
- Premium 9815



## Foot Prop Positioners

Maintains prescribed foot position. Facilitates offloading of the heel from the mattress surface.

- Premium Only
- Small 9825
- Small Flat-Sided 9826
- Large 9827



## Arm Cradles

Offers cushioning for arms, shoulders, elbows, and ulnar nerves. Design allows for varying arm positions.

- Standard 550236-CC
- Premium 9808



## Abduction Pillows

Assists in keeping patient's legs in proper position. Straps secure the pillow firmly. Optional comfort band available.

- Standard Medium 8015
- Premium Medium 9820

*Additional sizes and options available.  
See ordering information.*



## Arthroscopic Well Leg Holder

Provides support and comfort for legs, knees and ankles. Additional leg positioning products available. See ordering information.

- Premium Only 550223



## Roll Positioners

Tremendous versatility. Deep foam envelops the patient and offers comfort while managing pressure. Useful in a multitude of surgical and post-surgical procedures.



Foam	5" Dia	Half Roll
● <b>Standard</b>	8705	C3007
● <b>Premium</b>	9805	n/a

## Body Aligners

Promotes proper hip and spinal alignment while supporting patient repositioning as part of an effective pressure ulcer management program. Provides greater stability than pillows.



Small	W/O Sleeve	Cotton Cover
● <b>Standard</b>	8218	18071
● <b>Premium</b>	9811	n/a
Large	W/O Sleeve	Vinyl Cover
● <b>Standard</b>	8017	18070
● <b>Premium</b>	9812	n/a

## Positioning Wedges

Stabilizes the patient securely in several positions without slipping or sliding. Holds position better than rolled towels or pillows. Also used to support arm or leg when lying on side or at foot of bed to prevent sliding in supine or gatched positions.



Positioning Wedges	7"	10"	12"
<b>Wedge with Cotton Cover</b>	C7	C10	C12
<b>Wedge Vinyl Cover</b>	8040	8060	8050
<b>Wedge with Nylon Cover</b>	8073	8074	8072
Positioning Wedges	18"	12"	12"
<b>Wedge without Cover</b>	8020	8030	9810*

\*Premium density

## O.R. Table and Arm Board Pads

O.R. Table and Arm Board Pads help maintain hemodynamics during surgery. Use on top of an operating room table to enhance pressure redistribution.



Arm Board /Utility Pads		
● <b>Standard</b>	8616	16" x 8" x 1 1/2"
● <b>Premium</b>	550228-CC	20" x 8" x 2"
O.R. Table Pads		
● <b>Standard Only</b>	2074	74" x 20" x 4"
	2072	72" x 20" x 2"

## Neonate Crate

The soft Eggcrate® BioFoam pad, with sidewall construction and uniquely designed material, gently cradles the infant in hundreds of soft uniformly shaped “pillows” to create a comfortable and secure environment. The result— pressure redistribution and a gentle protection of the infant’s skin and tender areas.

- Pressure redistribution
- Sidewall construction helps position arms and legs close to the body
- Thermoplastic drape conforms to the pad and contours to wick away moisture
- May be used on top of existing mattress or as a mattress replacement
- Designed to allow placement of 15" BioFoam wedge to raise infant’s head and torso area to assist in reducing the risk of regurgitation after feeding



Neonate Crate	550001
Neonate Wedge	550002
Neonate Comforter	550003

## BioClinic Stretcher Mattresses

The BioClinic stretcher mattresses upgrade your facility transportation mattresses to therapeutic support surfaces at a price that offers a savings over standard mattress replacements.

### BioClinic stretcher benefits include:

- **Comfort** — select from high-resiliency foam that molds to the body contours or our premium viscoelastic foam that envelops and supports the body while redistributing pressure over the entire surface.
- **Choice** — choose premium or standard foam density
  - Premium: 3.3 density viscoelastic topper, 1.7 density premium performance base
  - Standard: 2.25 density high-resiliency foam topper, 1.7 density premium performance base
- **Custom fit** — all mattresses are custom-made to fit your fleet of stretchers.
  - Available in three-inch or four-inch depths
  - Maximum dimensions 30" W x 80" L in .5" increments
- **Cost effective solutions** — at a price less than the cost of an original equipment replacement mattress, you can increase your therapeutic benefit while reducing your patient’s risk for skin breakdown



### Head and Body Positioners - Head Cradle Positioners

Part No	Description	Foam	Size	Qty
8813	Pediatric	Standard	4" x 5" x 2"	48/cs
9813	Pediatric	Premium	4" x 5" x 2"	48/cs
8814	Child	Standard	7" x 6" x 3"	24/cs
9814	Child	Premium	7" x 6" x 3"	24/cs
8815	Adult	Standard	9" x 8" x 4"	24/cs
9817	Adult	Premium	9" x 8" x 4"	24/cs
9824	Adult, Slotted	Premium	9" x 8" x 4"	24/cs
9818	Adult, Large Frame	Premium	11" x 10" x 5"	24/cs
550269	Adult, Large Frame Visco	Premium	11" x 10" x 5"	20/cs

### Body Aligners

Part No	Description	Foam	Size	Qty
8218	W/O Sleeve - Small	Standard	7" x 7" x 17"	8/cs
18071	W/Cotton Cover - Small	Standard	7" x 7" x 17"	8/cs
9811	W/O Sleeve - Small	Premium	7" x 7" x 17"	8/cs
8017	W/O Sleeve - Large	Standard	9 7/16" x 8" x 23"	4/cs
18070	W/Vinyl Cover - Large	Standard	9 7/16" x 8" x 23"	4/cs
9812	W/O Sleeve - Large	Premium	9 7/16" x 8" x 23"	4/cs

### O.R. Foam Protectors - Convoluted Arm Board/Utility Pad

Part No	Description	Foam	Size	Qty
8616	Arm Board Pad	Standard	16" x 8" x 1"	12 pr/cs
550228-CC	Arm Board Pad	Premium	20" x 8" x 2"	12 pr/cs

### Positioning Wedges

Part No	Description	Foam	Size	Qty
8040	7" Wedge w/Vinyl Cover	Standard	24" x 24" x 7"	2/cs
8073	7" Wedge w/Nylon Cover	Standard	24" x 24" x 7"	2/cs
C7	7" Wedge w/Cotton Cover	Standard	24" x 24" x 7"	2/cs
8060	10" Wedge w/Vinyl Cover	Standard	24" x 24" x 10"	2/cs
8074	10" Wedge w/Nylon Cover	Standard	24" x 24" x 10"	2/cs
C10	10" Wedge w/Cotton Cover	Standard	24" x 24" x 10"	2/cs
8050	12" Wedge w/Vinyl Cover	Standard	24" x 24" x 12"	2/cs
8072	12" Wedge w/Nylon Cover	Standard	24" x 24" x 12"	2/cs
C12	12" Wedge w/Cotton Cover	Standard	24" x 24" x 12"	2/cs
8020	18" Wedge w/o Cover	Standard	18" x 25" x 18"	2/cs
8030	12" Wedge w/o Cover	Standard	30" x 21" x 12"	2/cs
9810	12" Wedge w/o Cover	Premium	20" x 21" x 12"	2/cs

### Round Foam Ring Positioners

Part No	Description	Foam	Size	Qty
8804	Round Ring Set	Standard	5" dia.	24pr/cs
8805	Round Ring Set	Premium	5" dia.	24pr/cs
8807	Round Ring Set	Standard	7" dia.	24pr/cs
8809	Round Ring Set	Standard	9" dia.	32/cs
550211-CC	Round Ring Set	Premium	9" dia.	32/cs
9809	Multi Set 3", 5", 7", 10"	Premium	1" dia.	32/cs

**O.R. Foam Protectors - O.R. Table Pad**

Part No	Description	Foam	Size	Qty
2074	4" Convuluted Foam Pad	Standard	72" x 20" x 4"	3/cs
2072	2" Convuluted Foam Pad	Standard	72" x 20" x 2"	12/cs

**Roll Positioners**

Part No	Description	Foam	Size	Qty
8705	Positioning Roll	Standard	5" dia. x 17"	9/cs
9805	Positioning Roll	Premium	5" dia. x 17"	9/cs
C3007	Half-Roll	Standard	5" x 17" x 7"	9/cs

**Foam Towel**

Part No	Description	Foam	Size	Qty
8212	Foam Towel	Standard	108" x 12"	12/cs

**Extremity Positioners - Ulnar Nerve Protectors**

Part No	Description	Foam	Size	Qty
4816	Ulnar Nerve Protector	Standard	6" x 6"	20 pr/cs
9816	Ulnar Nerve Protector	Premium	6" x 6"	20 pr/cs

**Laminectomy Arm Cradles**

Part No	Description	Foam	Size	Qty
550236-CC	24" x 5" x 3"	Standard	24" x 5" x 3"	12 pr/cs
9808	24" x 5" x 3"	Premium	24" x 5" x 3"	12 pr/cs

**Arthroscopic Knee Holder**

Part No	Description	Foam	Size	Qty
9829	Knee Holder	Premium	36" x 1" x 6"	8/cs

**Abduction Pillows**

Part No	Description	Foam	Size	Qty
8022	Contoured Small	Standard	12" x 18"	4/cs
8122	Contoured w/Comfort Band Small	Standard	12" x 18"	4/cs
9819	Contoured Small	Premium	12" x 18"	4/cs
9719	Contoured w/Comfort Band Small	Premium	12" x 18"	4/cs
8037H	Flat Sided Small	Standard	12" x 18"	4/cs
8015	Contoured Medium	Standard	14" x 22"	4/cs
8115	Contoured w/Comfort Band Medium	Standard	14" x 22"	4/cs
9820	Contoured Medium	Premium	14" x 22"	4/cs
9720	Contoured w/Comfort Band Medium	Premium	14" x 22"	4/cs
8035H	Flat Sided Medium	Standard	14" x 22"	4/cs
8023	Contoured Large	Standard	16" x 25"	4/cs
8123	Contoured w/Comfort Band Large	Standard	16" x 25"	4/cs
9821	Contoured Large	Premium	16" x 25"	4/cs
9721	Contoured w/Comfort Band Large	Premium	16" x 25"	4/cs
8036H	Flat Sided Large	Standard	16" x 25"	4/cs
8024	Contoured XL	Standard	20" x 22"	4/cs
8124	Contoured w/Comfort Band XL	Standard	20" x 22"	4/cs
9722	Contoured w/Comfort Band XL	Premium	20" x 22"	4/cs
8038H	Flat Sided XL	Standard	20" x 22"	4/cs

### Arthroscopic Well Leg Holder

Part No	Description	Foam	Size	Qty
550223	Well Leg Holder	Premium	20" x 4" x 10"	4/cs

### Leg Elevators

Part No	Description	Foam	Size	Qty
8016	Regular	Standard	11" x 15" x 8"	4/cs
8021	Long	Standard	11" x 31" x 14"	1/cs

### Foot Prop Positioners

Part No	Description	Foam	Size	Qty
8032	Convuluted Eggcrate	Standard	6" x 22" x 4"	4/cs
9825	Small	Premium	8" x 7" x 8"	6/cs
9826	Small, Flat-Sided	Premium	8" x 7" x 8"	6/cs
9827	Large	Premium	8" x 8" x 9"	4/cs

### Foot Positioners

Part No	Description	Foam	Size	Qty
8028	Foot Positioner	Standard	14" x 8" x 8"	4/cs
9828	Foot Positioner	Premium	14" x 8" x 8"	4/cs

### Heel/Ankle Protectors

Part No	Description	Foam	Size	Qty
4815	Heel/Ankle Protector	Standard	6" x 6"	20 pr/cs
9815	Heel/Ankle Protector	Premium	6" x 6"	20 pr/cs

### Neonate Crate

Part No	Description	Size	Qty
550001	Neonate Crate (Includes Eggcrate Pad, Comforter and protective drape)	12" x 25" x 2"	10/cs
550002	Neonate Wedge	11" x 12" x 3"	20/cs
550003	Neonate Comforter (Facing 75% polyester/25% acrylic; backing 100% olefin) Passes flammability testing/non-toxic	25" x 12"	20/cs

