



Joerns
Healthcare

Joerns® P.R.O. Matt Plus Non-Powered Convertible Air Foam Therapy System



transforming wound management



Mattress Features and Benefits

- Non-powered mattress replacement converts to alternating pressure therapy with addition of optional control unit
- Four-zone system (head, shoulder, torso and foot) inner air core reproduces the surface pressure profile of a therapeutic low air loss mattress
- P.R.O. Technology valve maintains optimal air pressures in each therapy zone for superior pressure redistribution regardless of patient weight and profile
- Patient is primarily on air mattress versus foam
- One-inch viscoelastic topper enhances pressure redistribution and patient comfort
- Surrounding firm perimeter provides stability during patient care and contributes to patient safety by eliminating gaps between frames and siderails
- Self-calibrating design needs no manual inflation
- Nylon therapy pad is standard and includes these unique qualities:
 - Waterproof
 - Low friction and low shear surface
- Stretch cover optional
- Maintenance-free with the exception of normal surface cleaning

Optional Control Unit Specifications

Feature	ES Control Unit	ES2 Control Unit	CBC Control Unit
Dimensions	13.5" (W) x 9"(H) x 8.5"(D)	7.5" (W) x 12.25"(H) x 5.5"(D)	11.18"(W) x 5.6"(H) x 3.5"(D)
Weight	11 lbs.	7 lbs.	4.2 lbs.
Comfort settings	8	5	10
Cycle times	5, 10, 15 minutes	5, 10, 15 minutes	10 minutes
Modes of Operation	Three - <i>Autofirm, Alternate and Therapy</i>	Three - <i>Autofirm, Alternate and Therapy</i>	Two - <i>Alternate and Static</i>
Lockout	Yes	Yes	No
Warranty	Two-years	One-year	One-year
Agency Approvals	UL Classified: UL 60601-1 CAN/CSA C22.2 No. 60601-1	UL Classified: UL 60601-1 CAN/CSA C22.2 No. 60601-1	UL Classified: UL 60601-1 CAN/CSA C22.2 No. 60601-1
Electrical Specifications	120V AC, 60 Hz, 0.6A	120V AC, 60 Hz, 0.6A	120V AC, 60 Hz, 0.6A

Clinical Indications

Pressure Redistribution: Stage I – IV Pressure Ulcers*; Rehabilitation; Burns; Neurology; Dermatology; Amputations

Pain Management: AIDS; Arthritis; Oncology

**In the non-powered mode, the P.R.O. Matt Plus is indicated for prevention through Stage II pressure ulcers. With the addition of an optional control unit and based on individual patient assessment, the P.R.O. Matt Plus may be used for the treatment of uncomplicated Stage III or IV pressure ulcers.*

Technical Specifications

Maximum Weight Capacity*: 500 lbs

Mattress Dimensions

- 36"(W) x 80"(L) x 7"(D)
- 36"(W) x 84"(L) x 7"(D)
- 42"(W) x 80"(L) x 7"(D)

Meets Flamability Standards

TB 117, TB 129, 16 CFR 1633 and BFD IX-1

Mattress Warranty

Five-years, non-prorated

Environmental Conditions

Operating Conditions

- Ambient Temperature: +10°C to +40°C
- Relative Humidity: 30% to 75% Non-Condensing

Storage and Shipping Conditions

- Ambient Temperature: -10°C to +40°C
- Relative Humidity: 10% to 100%

**Mattress weight capacity only; total weight must not exceed bed frame manufacturer's specified load capacity. Joerns reserves the right to change specifications without notice.*