



## Features and Benefits

### Control Unit Therapy Features

- Low air loss therapy addresses microclimate control and provides pressure redistribution
- Three modes of operation -Autofirm, Therapy and Pulsate
  - *Autofirm* provides maximum air inflation to assist both resident and caregiver during resident care and transfer
  - *Therapy* control settings with eight comfort adjust levels to maximize resident comfort and promote healing
  - *Pulsate* alternates between a decrease and increase of the inflation in all cells every 30 seconds
- Automatic panel lockout to avoid unintended adjustments
- Audible and visual power failure alarm



### Mattress Therapy Features

- Twenty therapy cells with eight-inch deep construction for superior pressure redistribution
- Therapy pad with unique qualities:
  - Fluid proof
  - Moisture vapor permeable
  - Low friction and shear surface
- Rapid CPR deflation with quick disconnect hose feature
- Lower 2-inch foam safety cell provides support in the event of a power failure or for patient transport

### Clinical Indications

**Pressure Relief:** Stage I – IV Pressure Ulcers; Moisture Management Required; Rehabilitation; Burns; Neurology; Dermatology; Amputations  
**Pain Management:** AIDS; Arthritis; Oncology

## Technical Specifications

**Maximum Weight Capacity\*:** 450 lbs

### Control Unit

Weight: 9.9 lbs

Dimensions: 7.5" (W) x 12.25" (H) x 5.5" (D)

### Mattress

Dimensions: 36" (W) x 80" (L) x 10" (D)

### HCPCS Code

E0277 (Group II)

### Warranty

Control Unit: 2 years

Mattress: 1 year

### Agency Approvals

UL Classified: UL 60601-1

CAN/CSA C22.2 601.1-M90

### Electrical Specifications

120V AC, 60 Hz

### Environmental Conditions

Operating Conditions

Ambient Temperature: 10°C to +40°C

Relative Humidity: Non condensing

Storage and Shipping Conditions

Ambient Temperature: 10°C to +40°C

Relative Humidity: Non-Condensing

\*Mattress weight capacity only; total weight must not exceed bed frame manufacturer's specified load capacity.  
 Joerns reserves the right to change specifications without notice.



www.joerns.com • email: info@joerns.com

### Post Acute, Acute, HomeCare

5001 Joerns Drive  
 Stevens Point WI 54481  
 (P) 800.826.0270  
 (F) 800.457.8827

### VA/Government

19748 Dearborn Street  
 Chatsworth, CA 91311  
 (P) 800.966.6662  
 (F) 800.232.9796

### Canadian Office

1000 Clarke Road, Suite 6  
 London, ON Canada N5V 3A9  
 (P) 866.546.1151  
 (F) 519.451.8662