

## Introducing the NEW! CareWide Extension Kit

### Preventing falls reduces injuries and saves your facility unnecessary expenses.

Resident falls from the bed have been documented as a prevalent occurrence in post acute facilities for several years. Additionally, the average healthcare cost of a fall injury is \$19,400 (Thizio, 2004). Recently, a new study suggests that facilities can reduce the fall risk by 14% by utilizing a 39-inch extension kit on a standard size bed frame (Perry, 2010).

Joerns Healthcare is pleased to introduce the new **CAREWIDE EXTENSION KIT** available for our Easy Care® 3, 5, 7 and 9 bed frames. The new CareWide extension kit will transform your Easy Care beds by adding four inches to the surface width to help with independent repositioning and provide additional comfort for the resident.



Easy Care 5 Shown

#### CareWide Extension Kit Features:

- Transforms the Easy Care 3, 5, 7, and 9 bed frames to 39" width from 35" width
- Easily attached with a minimum of tools to mattress support platform
- The CareWide Extension Kit accommodates all Joerns® panel styles and finishes
- Joerns assist devices work in tandem with the CareWide accessory (*includes the F14SC Assist Device, F028 Deluxe Assist Handle and the F025 Assist Handle*)
- Includes 39" wide panels and 39" mattress stops
- Compatible with Easy Care 3, 5, 7, and 9 bed frames currently in the field

Our NEW! CareWide Kit helps redefine your care environment while delivering safe and more comfortable care. Contact your Account Executive to order your CareWide Extension Kit today! Call 800.826.0270.