

Introducing Joerns BCX

Redefining the Bariatric Experience

The bed that harmonizes safe care with dignity and comfort



Joerns BCX
Bariatric Care



Designed for Bariatric Resident Safety and Security

Ensure uncompromised safety, dignity and comfort through advanced engineering that delivers best-in-class low bed height, a fully customizable patient environment, and ergonomic right-sized side rails

Engineered for Bariatric Care without Compromise

Maximize patient comfort with an intelligent bed designed for individualized positioning, featuring tool-less frame adjustments



Unique Platform Design

Enhance operational efficiencies with an intuitive interface and open platform architecture, designed to maximize asset utilization, lower cost of ownership and simplify maintenance

Call your representative or Joerns Customer Care at 800.826.0270 to schedule a demonstration now!

BCX Specifications

Overall Length: 84.25" - 92.25"

Overall Width: 41" - 53"

Weight Capacity: 750 lb.

Height Range: 8" - 32.5"

Caster Diameter:
4" AutoLock caster system

Mattress Size:
Width 36" - 48"
Length 80" - 88"

Mobility Height: 12" programmable

Side Assist: Tool-less Removal

Controls: Advanced Comfort Positioning