



## V.A.C. SIMPLICITY™ THERAPY SYSTEM



# V.A.C. Simplicity™ Negative Pressure Wound Therapy System

Wound healing at the touch of a button

### Easy to Operate

- Simple interface minimizes device training time
- Automatically delivers the most commonly prescribed NPWT setting of Continuous -125 mmHg
- SensaT.R.A.C.™ Technology continually monitors and regulates prescribed pressure at the wound site
- Quick Release canisters allow for efficient replacement
- Multi-orientational canister allows for patient mobility and product position flexibility
- 300mL canister minimizes canister changes
- Long battery life (14+ hours) with fast recharge time

### Patient Comfort

- Lightweight and portable
- Carrying case allows discreet delivery of therapy on-the-go
- Alarms help provide patient safety

### Specifications:

V.A.C. Simplicity™ Therapy Unit: P/N 414063  
Dimension: 7.6" (19.3cm) W x 6" (12.2cm) H x 6.4" (6.4cm) D  
Weight: 2.3 lbs. (1.04 kg)  
Canister volume: 300 ml  
Battery Life: 14+ hours (recharge time <6 hours)  
Electrical Charger/External Supply Input:  
100-240 VAC .72 A @ 115 VAC 50/60 Hz  
Charger/External Supply Output: 12 V, 3.3 A  
Maximum Electrical Leakage: <100 Microamps  
IEC Classification: Type B Applied Part, Class II, IPX0  
Specifications are subject to change without notice.

### Therapy Options:

Negative Pressure: -125 mmHg  
Mode of Operation: Continuous  
Intensity: High

### Storage Conditions:

Temperature Range ..... -4°F (-20°C) to 140°F (60°C)  
Relative Humidity Range .... 0-95% non-condensing

### Operating Conditions:

Temperature Range ..... 41°F (5°C) to 104°F (40°C)  
Relative Humidity Range .... 0-95% non-condensing



**NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for KCI products and therapies. Please consult product instructions for use prior to application. Rx only.**

©2012 KCI Licensing, Inc. All rights reserved. Joerns is a trademark of Joerns RecoverCare. All other trademarks designated herein are proprietary to KCI Licensing Inc., its affiliates and/or licensors. DSL#15-0245.US.JRC (Rev. 4/15) 6240001 RevA • 15-2201

