

## Expandable Fall Prevention Bed

Minimizing risk of patient falls can save a facility thousands or even millions in costs and improve the standard of care. The RC Ultra Hi-Lo 850 lowers to 9½" and has total lock and steering lock casters to allow for easy bed movement and facilitate patient transport. The simple to use pull pin system allows caregivers to quickly and easily expand the bed frame's width to accommodate various patient girths and emergency transport through 40" door frames. Unlike most hi-lo beds, the RC Ultra Hi-Lo 850 can roll at any height making it safer for the caregivers to move throughout the facility.



## A fall prevention intervention meeting the highest standards

**WEIGHT MONITORING**

Elevate safe patient handling with weigh-in-bed capability to lower the risk of falls and skin shear

**BED EXIT WITH CONNECTIVITY**

Enhance patient monitoring and response with nurse call connectivity and bed exit alarm

**LEADING LOW HEIGHT**

Reduce falls and injuries with 9½" low height and allow for more patients to actively ambulate

## Standard Features

- Descends to a minimum deck height of 9½"
- 5 functions including Trendelenburg and Comfort Chair
- Integrated foot end staff control with patient lock-out
- Bed deck expands to 39", 42" or 48"
- Durable, 10 button hand pendant
- Total lock and steering lock casters for easy movement at any height
- Optional integrated scale for one-button weighing in any position
- Optional integrated bed exit and nurse call connectivity
- Optional trapeze

## Frame Specifications

<b>Bed Deck Length:</b>	84"
	(96" with optional extender)
<b>Overall Frame Size:</b>	36" to 48" W to 89" L
<b>Mattress Size:</b>	Minimum 36"
	Maximum 48"
<b>Bed Deck Height:</b>	Minimum 9½"
	Maximum 29"
<b>Head/Foot Panel Width:</b>	38"
<b>Safe Working Load:</b>	850 lbs
<b>Power Requirements:</b>	115 VAC, 60 Hz, 10.5 AMP
<b>Warranty:</b>	Structural 5 years
	Electrical System 1 year

## RC Ultra Safety Suite

### Includes:

- RC Ultra Hi-Lo 850
- Pressure redistribution mattress with Safety Perimeter
- Bed Alarm for patient exit notifications
- Fall pad for added piece of mind



**UPGRADE THE SUITE** with a lift or CareChair for safe patient handling