

RecoverAir 2000 Elite

*Low Air Loss and
Alternating Pressure Overlay*



Dual therapy features for advanced wound treatment

Powered low air loss provides excellent pressure redistributing therapy and creates air flow under your patients, protecting against maceration. The low friction, vapor permeable comforter protects from shearing and provides a quilted “wicking” material underneath to help keep your patients dry. The RecoverAir 2000 Elite overlay utilizes a unique 3-1 air cell inflation/deflation cycle for alternating therapy, assuring optimum pressure redistribution and clinical outcomes. Patient comfort settings allow adjustments so your patients remain as comfortable as possible. The RecoverAir 2000 Elite overlay is complete care for your patients.

RecoverCare[®]
Equipment. Solutions. Results!

RecoverAir 2000 Elite *Low Air Loss and Alternating Pressure Overlay*

Features

- Provides effective low air loss to assist in moisture control
- Selectable 3 - 1 air cell alternating pressure therapy
- Maximum inflate for easy patient positioning
- Low friction, anti shear cover
- CPR quick release for rapid deflation
- Whisper quiet operation
- Latex free

Clinical Considerations

- Prevention and treatment of Stage I - IV pressure ulcers, pain management, immobility
- Contra-indicated for patients with unstable spine conditions



Specifications

- Dimensions
 - Pump: 14" x 5" x 9" (13 lbs)
 - Overlay: 80" x 35.5" x 3.5"
- Weight limit: 500 lbs
- Electrical rating: 120V~, 60Hz
- Agency Approvals: UL; BSI; TUV

For Ordering Information

PHONE 888-750-7828 • FAX 866-750-7828

www.recovercare.com

1920 Stanley Gault Parkway, Suite 100, Louisville, KY 40223

RecoverCare[®]
Equipment. Solutions. Results!