

Joerns ACX Active Care Bed: Expanding Recovery Care Solutions

The Joerns ACX bed is designed to expand to fit more patient care plans. Width has been expanded. Capacity has been expanded. Technology has been expanded. Critical care can be expanded.

The ACX Active Care bed offers a complete care solution in an expandable low-bed to enhance patient mobility, expand use to more patients and efficiently monitor patients.



WEIGHT MONITORING




Elevate safe patient handling with weigh-in-bed capability to lower the risk of falls and skin shear

BED EXIT WITH CONNECTIVITY



Enhance patient monitoring and response with nurse call connectivity and bed exit alarm

INTEGRATED WIDTH IN CRITICAL CARE



Expandable width for patient comfort and repositioning while enabling emergency transport

Expand the dimension of care

Key Features

Advanced care

- 600lb capacity of high utilization
- Expandable to 48" width and 88" length for comfort and fit
- Easy-view side rails allow for improved patient visibility with convenient hand holds for active repositioning
- Advanced Comfort Chair and contour positioning for reduced skin shear and improved comfort
- Regression articulated back section to reduce pressure when raising head of bed

Safety

- Low-bed height of 10" for early patient mobility and falls management
- In-rail controls for attentive bedside care
- Electric and manual CPR release
- Safety-stop technology for obstruction-free operation
- Foot end and in-rail digital back angle display for confirmed patient positioning
- Under bed light for safe patient mobility
- Optional scale for weight monitoring
- Optional nurse call and 3-zone bed exit alarm for earlier fall detection

Mobility

- Center lock and steer casters for efficient mobility
- Integrated therapy pump holder
- Optional on-board battery for enhanced transport capability
- Optional integrated hospital-grade outlet for transport convenience

Expand the dimension of care while reducing falls risk, enabling safe patient handling and effectively monitor patient safety



Specifications

Overall Length	88.5"
Overall Width	42"
Weight Capacity	600 lb
Weight	Basic: 418 lb Premium: 466 lb
High Height	31"
Low Height	10"

Mattress Size	36" - 48" Wide 80" - 88" Long
Positioning	70° Head angle 35° Knee angle 14° Trend / Reverse Trend
Caster Diameter	5"
Power Requirements	115 VAC, 60 Hz, 4 AMP