

# PrevaMatt® Series

Therapeutic Support Surfaces



# Superior Technology = Superior Performance

The Joerns Healthcare PrevaMatt® series has evolved into a portfolio of foam mattresses designed for comfort and safety all while providing pressure injury prevention. The focused portfolio gives healthcare Providers a range of options to improve outcomes, reduce risks to both the caregiver and patient, and meet financial objectives.

- **IMPROVE CLINICAL OUTCOMES**  
Excellent pressure redistribution aids in wound prevention
- **REDUCE RISK OF INJURY**  
Standard safety features reduce fall risk and promote mobility
- **REDUCE FINANCIAL BURDEN**  
Broad and diverse portfolio designed for quality and durability



GREENER & MORE  
SUSTAINABLE PORTFOLIO

## Environmentally Focused

- Compatible with recovery and recycling technologies and local requirements
- CertiPUR-US® certified foams improve the health, environmental and safety experience of patients

## Sustainability

- Design and content enhance useful life improving consumption and sustainability
- Elimination of unnecessary chemicals results in purer material and can lead to improved durability



## MANIPULATION OF FOAM PROPERTIES



### Open cell memory foam technology

- Allows greater air flow resulting in more breathability and a cooler, drier patient experience
- Materials less affected by fatigue – better height and support retention with faster recovery time

### Surface Modification and Layering

- Alterations in pattern, size of pattern, spacing between patterns and depth of cut to protect high risk areas
- Zoned pressure redistribution with improved air circulation and maximized comfort



## TEMPUR® Option



Partnership with Tempur-Pedic® to provide genuine TEMPUR® material on the PrevaMatt Alleviate

- Reduces wound incidence rates
- Treats up to Stage 3 uncomplicated pressure injuries



# SAFETY AND SECURITY

## Firm Perimeter



















- Provides edge of bed stability
- Aids ingress and egress

## Optional Raised Perimeter

- Aids in fall prevention
- Reduces risk of entrapment

When used as part of a Joerns bed system, our PrevaMatt mattresses are guaranteed to meet the guidelines set forth in the Hospital Bed System Dimensional and Assessment Guidance to Reduce Entrapment developed by the Food and Drug Administration (FDA) and the Hospital Bed Safety Workgroup (HBSW).

## PrevaMatt® Series of Prevention Surfaces

	BASIC	DEFEND	CONSOLE	ALLEVIATE
SUSTAINABLE MATERIAL				
SURFACE MODIFICATION + LAYERING				
NATURAL INGREDIENTS				
COVER MATERIAL (NYLON, 2-WAY, 4-WAY STRETCH)				
OPEN CELL VISCO				
TEMPUR® OPTION				
CLINICAL INDICATIONS	Prevention	Prevention to Stage 1	Prevention to Stage 1	Prevention to Stage 2 <sup>†</sup>

<sup>†</sup>TEMPUR® option is suitable for treatment of up to Stage 3 uncomplicated pressure injuries

<sup>\*\*</sup>TEMPUR is a registered trademark of Tempur-Pedic®