

Expandable Bariatric Bed Frame

Is an increase in morbidly obese patients causing your budget to explode? The RC 750 is a 750 lb capacity bed frame with three function electric control. The simple to use pull pin system allows caregivers to easily expand the bed frame's width from 39" to 54" to accommodate a variety of body contours and emergency transport through 40" door frames. The RC 750's economical price makes it an excellent bariatric bed solution for the alternate care market.



Multi-functional Bariatric Capacity in an Affordable Package

MAXIMUM COMFORT ----54"--- 80" 96" Greatest range of comfort with 54" Width and 96" Length capability







Standard Features

- · Split frame design for easy delivery and storage
- Expandable bed deck permits transport through doorways as narrow as 40"
- Accommodates mattresses up to 54" wide
- 3 position, wing type, padded side rails
- 3 functions including hi-low, head (Fowler) and knee gatch
- · IV pole and Foley bag holder
- 750 pound weight capacity
- · Optional trapeze

Frame Specifications

| • | | |
|--|---------------------------|----------------|
| Bed Deck Length: | 80" (96" with optic | onal extender) |
| Overall Frame Size: | 39.9" to 56.5" W to 88" L | |
| Mattress Size: | Minimum | 37" |
| | Maximum | 54" |
| Bed Deck Height: | Minimum | 17" |
| | Maximum | 30" |
| Head/Foot Panel Width: | | 38" |
| Safe Working Load: | | 750 lbs |
| Power Requirements: 115 VAC, 60 Hz, 10.5 AMP | | |
| Warranty: | Structural | 5 years |

RC 750 Suite Package

Includes:

- RC 750
- Pressure redistribution mattress with Safety Perimeter
- Fall pad for added piece of mind

UPGRADE THE SUITE with a lift or CareChair for safe patient handling





Electrical System

1 year