



Active Transfer Device for Use in All Settings

Easy to use and non-powered, Hoyer Up is a manual stand-aid and seated transfer device suitable for sling-free transfers with caregiver assistance.

Hoyer Up is designed for users who have some weight-bearing ability and can participate in the transfer promoting patient rehabilitation and functional independence. The Up simplifies care in all settings and helps maintain more activities of daily living.



UNIQUE MODULAR DESIGN




Simplifies inventory and transport with an efficient flat-pack design that requires no tools for assembly

SWING SEAT




Secure padded swing seat surface delivers comfort and stability. Seats fold flat against the side arms when not in use.

COMFORTABLE KNEE PAD



Enhances comfort and accommodates a larger range of patient sizes with a contoured, height adjustable knee pad.

ADJUSTABLE BASE



Increases patient care accessibility with a simple and robust mechanical foot pedal providing open/close functionality.



SAFE WORKING LOAD: 440 LB / 200 KG

SPECIFICATION

Overall Length	37.4"	950 mm
Overall width	24.8"	630 mm
Overall Height	43.3"	1100 mm
Legs Open - External Width	33.3"	845 mm
Legs Open - Internal Width	28.5"	725 mm
Legs Closed - External Width	23.6"	600 mm
Legs Closed - Internal Width	18.9"	480 mm
Leg Height	4.3"	110 mm
Casters	3" Front	4" Rear

WEIGHTS

Total*	72.8 lb	33 kg
--------	---------	-------

*Base 43 lbs, seat 21 lbs, mast 8.8 lbs

Specifications may vary due to manufacturing tolerances.