

Laminate Bed Panels

We have panel style and finish choices to complement your décor and furnishings. When selecting bed panels, special consideration should be given to compatibility, durability, safety and resident entrapment concerns. Joerns Healthcare recommends the use of Joerns bed panels with Joerns beds when configuring a bed system.

Panel Options¹

Classical



DIMENSIONS²
36" Standard: 36"W x 0.75"D x 18"H
42" Wide: 42"W x 0.75"D x 18"H

Traditional



DIMENSIONS²
36" Standard: 36"W x 0.75"D x 18.25"H
42" Wide: 42"W x 0.75"D x 18.25"H

Optional Staff Control

UltraCare[®] and EasyCare[®]




WeCare[®]



Color Options



Natural Maple with Sand Edge



Quartered Oak with Sand Edge



Spiced Cherry with Brown Edge



Bayou Walnut with Brown Edge



Cordovan with Brown Edge

Universal Joerns Bed Panels³

Joerns' panel selection is now simpler than ever. Our universal bed panels are designed to fit all current and legacy Joerns bed models, making it easy to refresh your environment today while staying aligned with future upgrades.

Each non-staff control panel is fully interchangeable between the head and foot of the bed, reducing complexity in selection, inventory, and replacement.

One panel. Any bed. Maximum flexibility.

1. Finishes and hardware are subject to change over time. We reserve the right to make changes consistent with design trends or for reasons beyond our control. Finish options have been reproduced as faithfully as processes allow.
2. Dimensions rounded up to the ¼ inch; Measurements taken from outermost points.
3. Some restrictions apply. Contact Joerns Healthcare for your bed panel / model compatibility.
4. Non-Joerns panel kits available upon request.